

# FREEWHEEL

ALL TERRAIN WHEELCHAIR ATTACHMENT



The FreeWheel is a single large pneumatic castor that after attaching to your foot-plate or bottom rails, allows you to easily cross rough ground. It lifts your castors over the ground and rides on the large wheel to provide a comfortable and stable ride, easily riding over gravel, grass and dirt. Fit some mountain bike tires to your wheelchair and you could easily cross sand and snow.

This system is adjustable to suit most wheelchairs. It can be attached quickly (within seconds) to your wheelchair, while an

easy to use clamp keeps it secure.

Gerry is shown here in the gravel pit on the Melrose testing range although rough ground isn't the only place these excel on. Once on the road or the footpath the ability to tighten the pivot allows the wheelchair to track straight, making it a good distance trainer.

# FREEWHEEL

ALL TERRAIN WHEELCHAIR ATTACHMENT

## STANDARD SPECIFICATIONS

Designed for use on most fixed or folding frame wheelchairs

Large 12.5" Pneumatic Tyre rolls easily over broken ground

Black Satin Powder-coat finish

Lightweight, it only weighs 2.35 kg

Easy to attach or remove in seconds from your wheelchair

Can be mounted on the back of your wheelchair when not in use with a fitting that is attached to the rear of your wheelchair back

Adjustable Castor return to center (tracking)

Adjustable Castor centering position (offset for cambered footpaths or roads)

Adjustable lift height (how high it raises your Front Castors)

1 Year Warranty

