

CHILD'S ADJUSTABLE FIXED FRAME WHEELCHAIR

HAWK GROWABLE



The Melrose Hawk Grow-able sits alongside its big brother as our main lightweight adjustable position wheelchair. Being able to grow in width and length allows for a longer period before having to move to a larger ride.

Adjustable position provides the ability to fine tune the seating

and balance point so the child can get a more dynamic ride as their confidence increases.

The chair is custom fabricated with aircraft grade aluminium and robust high quality components, that should out last even the most active child.



HAWK GROWABLE

CHILD'S ADJUSTABLE FIXED FRAME WHEELCHAIR



STANDARD SPECIFICATIONS

Aluminium Frame 90420G

Seat Width: 14" Grow-able

Seat Depth: 13" - 16"

Adjustable Back Height and Angle

Folding Back

Adjustable Foot-plate Height

Adjustable Seat Height

Adjustable Balance Point

Castors: Roller-blade
4", 5", 6" Solid
4", 5", 6" Wide Solid

Quick release Spoked Wheels on Alloy Rims

Tyres: High Pressure Tread or Smooth

Brakes: Push to Lock Brakes
Alloy Scissor Brakes

Black Fabric Sling and Back with Leather Trim

Leg Strap

Five Year Frame Warranty

OPTIONS AND ACCESSORIES

Castors: 4", 5", 6" Solid Alloy
4", 5", 6" Soft Roll
6", 8" Pneumatic

Frog Leg Forks

Quad Release Axles

Wheels: Spoked 20", 22", 24"
BMX Wheels
Xcore Wheels
Spinergy Wheels

Push-rims: Rubber Coated - Red or Black
Vinyl
Foam
CarboLife
Natural Fit

Push handles: Removable
Fold Down

Brakes: Composite Scissor Brakes
Ds-Locks Brakes
Hill Climber Brakes
Brake Extension Levers

Anti-Tips with Plastic Wheels or Roller-blades

Spoke Guards: Clear or Black

Armrest: Desk or Full Length Pad
T Adjustable Height, with/without Skirt Guards

Adjustable Back Upholstery

Skirt Guards: Fabric Skirt Guards
Fabric covered ABS Skirt Guards
ABS Drop In Skirt Guards
Alloy or Carbon Fibre Skirt Guards

Under Chair Bag

Fabric or ABS Foot-plate Cover

Cushion, Impact Guards, Castor Covers

Interchangeable Camber Bar 0°, 3°, 6°